
Recombinant Human IGF2BP2 Protein

Catalog Number: 230-00573

Data Sheet / Certificate of Analysis (CoA)

Last Revised: 11/26/2024

Source

Species	Human
Accession Number	Q9Y6M1
Gene Symbols	IGF2BP2
Synonyms	Insulin-like growth factor 2 mRNA-binding protein 2, IGF-II mRNA-binding protein 2
Expressed Region	Met2-Thr220


Preparation

Expression System	E. coli
Tag	N-terminal 6x histidine tag
Purity	> 95%
Molecular Weight	26.9 kDa
Endotoxin Level	Not determined

Specifications

Format	Lyophilized powder
Formulation	Lyophilized from a 0.2 µm filtered solution in 40 mM Tris-HCl (pH 8.0), 50 mM NaCl.
Reconstitution	Briefly spin the vial and bring the contents to the bottom prior to opening. It is recommended to reconstitute at 0.5 - 1.0 mg/mL with sterile deionized water or 1x PBS.

SDS-PAGE Image



Shipping

The product is shipped with ice packs. Upon arrival, immediately store it at the temperature recommended below.

Storage/Stability

- Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
- Upon arrival, the lyophilized protein may be stored for 2 weeks at 4°C.
- For long term storage, it is recommended to store desiccated below -20°C in a manual defrost freezer.
- Generally, the shelf life is up to 12 months from date of receipt at -20°C or -80°C under sterile conditions.
- Following reconstitution, the protein may be stored for 2 weeks under sterile conditions at -20°C. For long term storage, it is recommended to make appropriate aliquots at -20°C or -80°C.