

---

## Recombinant Human AGR2 Protein

**Catalog Number:** 230-00596

**Data Sheet / Certificate of Analysis (CoA)**

Last Revised: 11/25/2024

### Source

---

<b>Species</b>	Human
<b>Accession Number</b>	O95994
<b>Gene Symbols</b>	AGR2
<b>Synonyms</b>	Anterior Gradient 2, HAG-2
<b>Expressed Region</b>	Arg21-Leu175

### Preparation

---


<b>Expression System</b>	E. coli
<b>Tag</b>	N-terminal 6x histidine tag
<b>Purity</b>	> 95%
<b>Molecular Weight</b>	20 kDa
<b>Endotoxin Level</b>	Not determined

### Specifications

---

<b>Format</b>	Lyophilized powder
<b>Formulation</b>	Lyophilized from a 0.2 µm filtered solution in PBS.
<b>Reconstitution</b>	Briefly spin the vial and bring the contents to the bottom prior to opening. It is recommended to reconstitute at 0.5 - 1.0 mg/mL with sterile deionized water or 1x PBS.

### SDS-PAGE Image



### Shipping

The product is shipped with ice packs. Upon arrival, immediately store it at the temperature recommended below.

### Storage/Stability

- Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
- Upon arrival, the lyophilized protein may be stored for 2 weeks at 4°C.
- For long term storage, it is recommended to store desiccated below -20°C in a manual defrost freezer.
- Generally, the shelf life is up to 12 months from date of receipt at -20°C or -80°C under sterile conditions.
- Following reconstitution, the protein may be stored for 2 weeks under sterile conditions at -20°C. For long term storage, it is recommended to make appropriate aliquots at -20°C or -80°C.