

Catalog Number: **230-00607**

Recombinant Human REG-1-alpha**Source**

Species	Human
Accession Number	P05451
Gene Symbol	REG1A
Gene ID	5967
Expressed Region	Gln23-Asn166
Synonyms	Lithostathine-1-alpha, ICRF, REG-1-alpha,


Preparation

Expression System	E.Coli
Tag	N-terminal 6x histidine tag
Purification	His-tag affinity purification by immobilized metal ion affinity chromatography (IMAC)
Purity	> 95%
Purity Determined By	SDS-PAGE under reducing conditions and visualized by Coomassie blue staining
Molecular Weight	18.6

Protein Specifications

Format	Liquid
Formulation	Supplied as a 0.2 um filtered solution in 4 M urea, 50 mM Tris-HCl (pH 8.0).
Concentration	Lot specific (see the label on the vial)
Preservative	None
Endotoxin Level	Not determined

SDS-PAGE Image



Shipping

Ice packs

Storage/Stability

Upon arrival, the protein may be stored for 2 weeks at 4 °C. For long term storage, it is recommended to store at -20 °C or -80 °C in appropriate aliquots. Avoid repeated freeze-thaw cycles.

This product is furnished for **LABORATORY RESEARCH USE ONLY**.

Not for diagnostic or therapeutic use.