
Recombinant Human PD-1 Ligand 2/PD-L2/CD273 Protein

Catalog Number: 230-01168

Data Sheet / Certificate of Analysis (CoA)

Last Revised: 11/12/2024

Source

Species	Human
Accession Number	Q9BQ51
Gene Symbols	PDCD1LG2
Synonyms	PD-1 ligand 2, PD-L2, PDCD1 ligand 2, Programmed death ligand 2, Butyrophilin B7-DC, B7-DC, CD273.
Expressed Region	Leu20-Thr220


Preparation

Expression System	E. coli
Tag	N-terminal His-tag
Purity	>90%
Molecular Weight	25 kDa

Specifications

Format	Liquid
Formulation	Supplied as a 0.2 µm filtered solution in 15 mM Na ₂ HPO ₄ (pH 7.4), 5 mM NaH ₂ PO ₄ , 0.5 M NaCl, 400 mM imidazole, 6 M urea.
Concentration	Refer to the label on the vial for the specific concentration, determined by Bio-Rad Bradford protein dye reagent concentrate.

SDS-PAGE Image



Shipping

The product is shipped with ice packs. Upon arrival, immediately store it at the temperature recommended below.

Storage/Stability

- Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
- Upon arrival, the liquid protein may be stored for 2 weeks at 4°C.
- For long term storage, it is recommended to store at -20°C or -80°C in appropriate aliquots.
- Generally, the shelf life is up to 6 months from the date of receipt at -20°C or -80°C under sterile conditions.