

Product Identification

Product Name Recombinant Mpox Virus IFN-alpha/beta receptor-like Secreted Glycoprotein, B16R
Catalog Number 230-01193

Source

Species Mpox Virus (MPXV)
Accession Number Q8V4R2
Gene Symbols B16R
Protein Name / Synonyms IFN-alpha/beta-receptor-like secreted glycoprotein, B16R
Expressed Region His21-Glu352


Preparation

Expression System Escherichia coli (E.coli)
Tag N-terminal histidine tag
Purification His-tag affinity purification by immobilized metal ion affinity chromatography (IMAC)
Purity >95%
Purity Determined By SDS-PAGE under reducing conditions and visualized by Coomassie blue staining
Molecular Weight 37 kDa

Protein Specifications

Format Liquid
Formulation Filtered solution in 50 mM Tris-HCl, 500 mM NaCl, 600 mM Arginine, 10% glycerol, pH 8
Concentration (lot specific) Lot specific (see the label on the vial), determined by Bradford protein assay.
Applications WB, ELISA

Images



Shipping Type

Blue ice

Storage

-20°C

This product is for research use only.

ISO 13485:2016 Certified

3607 Parkway Lane, Suite 100

Peachtree Corners, GA 30092

770-729-2992 | 1-888-494-8555

www.RayBiotech.com | info@raybiotech.com

©2024 RayBiotech Life, Inc.

