

Product Identification

Product Name Recombinant Mpox Virus Secreted IL-18 Binding Protein D6L
Catalog Number 230-01194

Source

Species Mpox Virus (MPXV)
Accession Number Q8V570
Gene Symbols D6L
Protein Name / Synonyms Secreted IL-18 binding protein, D6L
Expressed Region Met1-Lys126



Preparation

Expression System Escherichia coli (E.coli)
Tag N-terminal histidine tag
Purification His-tag affinity purification by immobilized metal ion affinity chromatography (IMAC)
Purity >95%
Purity Determined By SDS-PAGE under reducing conditions and visualized by Coomassie blue staining
Molecular Weight 14 kDa

Protein Specifications

Format Liquid
Formulation 50 mM Tris-HCl, 0.5 M NaCl, 0.4 M arginine, 1% mannitol, 5% trehalose, pH 8
Concentration (lot specific) Lot specific (see the label on the vial), determined by Bradford protein assay.
Applications WB, ELISA

Images



Shipping Type

Blue ice

Storage

-20°C

This product is for research use only.

ISO 13485:2016 Certified

3607 Parkway Lane, Suite 100

Peachtree Corners, GA 30092

770-729-2992 | 1-888-494-8555

www.RayBiotech.com | info@raybiotech.com

©2024 RayBiotech Life, Inc.

