
Recombinant Human CNTF Protein

Catalog Number: 230-01211

Data Sheet / Certificate of Analysis (CoA)

Last Revised: 11/12/2024

Source

Species	Human
Accession Number	P26441
Gene Symbols	CNTF
Synonyms	CNTF
Expressed Region	Ala2-Met200


Preparation

Expression System	E. coli
Tag	N-terminal His-tag
Purity	>90%
Molecular Weight	30 kDa

Specifications

Format	Liquid
Formulation	Supplied as a 0.2 μ m filtered solution in 40 mM Tris-HCl, 50 mM NaCl, 200 mM Arginine, 0.01% Tween 80.
Concentration	Refer to the label on the vial for the specific concentration, determined by Bio-Rad Bradford protein dye reagent concentrate.

SDS-PAGE Image



Shipping

The product is shipped with ice packs. Upon arrival, immediately store it at the temperature recommended below.

Storage/Stability

- Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
- Upon arrival, the liquid protein may be stored for 2 weeks at 4°C.
- For long term storage, it is recommended to store at -20°C or -80°C in appropriate aliquots.
- Generally, the shelf life is up to 6 months from the date of receipt at -20°C or -80°C under sterile conditions.