
Recombinant Mouse Tryptase Beta-2

Catalog Number: 230-01251

Data Sheet / Certificate of Analysis (CoA)

Last Revised: 11/12/2024

Source

Species	Mouse
Accession Number	P21845
Gene Symbols	Tpsb2
Synonyms	EC:3.4.21.59, Tryptase-2, Mast cell protease 6 (mMCP-6).
Expressed Region	Ile32-Ser276


Preparation

Expression System	E. coli
Tag	N-terminal His-tag
Purity	>90%
Molecular Weight	28 kDa

Specifications

Format	Liquid
Formulation	Supplied as a 0.2 µm filtered solution in 0.1 M NaHCO ₃ , 0.5 M NaCl (pH 8.0).
Concentration	Refer to the label on the vial for the specific concentration, determined by Bio-Rad Bradford protein dye reagent concentrate.

SDS-PAGE Image



Shipping

The product is shipped with ice packs. Upon arrival, immediately store it at the temperature recommended below.

Storage/Stability

- Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
- Upon arrival, the liquid protein may be stored for 2 weeks at 4°C.
- For long term storage, it is recommended to store at -20°C or -80°C in appropriate aliquots.
- Generally, the shelf life is up to 6 months from the date of receipt at -20°C or -80°C under sterile conditions.