
Recombinant Mouse Melanotransferrin (CD228) Protein

Catalog Number: 230-01258

Data Sheet / Certificate of Analysis (CoA)

Last Revised: 11/12/2024

Source

Species	Mouse
Accession Number	Q9R0R1
Gene Symbols	Melft
Synonyms	Membrane-bound transferrin-like protein p97 (MTf), CD228.
Expressed Region	Val20-Cys709


Preparation

Expression System	E. coli
Tag	N-terminal His-tag
Purity	>90%
Molecular Weight	76 kDa

Specifications

Format	Liquid
Formulation	Supplied as a 0.2 µm filtered solution in 20 mM Tris-HCl (pH 8.0), 50 mM NaCl.
Concentration	Refer to the label on the vial for the specific concentration, determined by Bio-Rad Bradford protein dye reagent concentrate.

SDS-PAGE Image



Shipping

The product is shipped with ice packs. Upon arrival, immediately store it at the temperature recommended below.

Storage/Stability

- Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
- Upon arrival, the liquid protein may be stored for 2 weeks at 4°C.
- For long term storage, it is recommended to store at -20°C or -80°C in appropriate aliquots.
- Generally, the shelf life is up to 6 months from the date of receipt at -20°C or -80°C under sterile conditions.