
Recombinant Human IL-12B Protein, HEK293 Lysate

Catalog Number: 230-10126-100

Data Sheet / Certificate of Analysis (CoA)

Last Revised: 11/26/2024

Source

Species	Human
Accession Number	P29460
Gene Symbols	IL12B
Synonyms	Interleukin 12B, NKSF2, Interleukin 12B (Natural Killer Cell Stimulatory Factor 2, Cytotoxic Lymphocyte Maturation Factor 2, P40), Cytotoxic Lymphocyte Maturation Factor 40 KDa Subunit, NK Cell Stimulatory Factor Chain 2, Interleukin-12 Beta Chain, IL-12 Subunit P40, CLMF P40, IL-12B, Natural Killer Cell Stimulatory Factor, 40 KD Subunit, Natural Killer Cell Stimulatory Factor-2, Cytotoxic Lymphocyte Maturation Factor 2, Natural Killer Cell Stimulatory Factor, Interleukin-12 Subunit Beta, Interleukin 12, P40, IL12, Subunit P40, Interleukin 12, 40 KD Subunit, Subunit P40, CLMF2, IMD28, IMD29, IL12, CLMF, NKSF, P40
Expressed Region	Ile23-Ser328

Preparation

Expression System	HEK293 cells
Tag	N-terminal 6x histidine tag
Purity	Unpurified cell lysate. HEK293 cells transfected with expression vectors harboring target gene were harvested and washed with PBS twice. The cell pastes were re-suspended with ice-cold PBS containing mammalian cell protease inhibitor cocktail and further lysed with freeze-thaw cycles. After clarifying with 20,000 g centrifugation at 4°C for 30 min, the lysate was aliquoted, lyophilized, and stored at -80°C immediately. Protein concentration was determined by BCA kit (Thermo Scientific, Inc.) using BSA as protein standard. The gene overexpression in lysate was confirmed by Western blotting using anti-His tag antibody and/or target-specific antibodies and the lysate derived from HEK293 cells transfected with the empty expression vector was used as a negative control.
Molecular Weight	Recombinant protein has a calculated molecular mass of 37. The actual molecular weight may increase slightly due to the potential post-modifications (PTMs).

Endotoxin Level Not determined

Specifications

Format Lyophilized powder

Formulation Lyophilized from a 0.2 μm filtered solution in PBS containing mammalian cell protease inhibitor cocktail

Concentration Determined by BCA protein assay kit (Thermo Scientific)

Reconstitution Briefly spin the vial and bring the contents to the bottom prior to opening. It is recommended to reconstitute at 0.5 - 1.0 mg/mL with sterile deionized water or 1x PBS.

Shipping

The product is shipped with ice packs. Upon arrival, immediately store it at the temperature recommended below.

Storage/Stability

- Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
- Upon arrival, the lyophilized protein may be stored for 2 weeks at 4°C.
- For long term storage, it is recommended to store desiccated below -20°C in a manual defrost freezer.
- Generally, the shelf life is up to 12 months from date of receipt at -20°C or -80°C under sterile conditions.
- Following reconstitution, the protein may be stored for 2 weeks under sterile conditions at -20°C. For long term storage, it is recommended to make appropriate aliquots at -20°C or -80°C.