

---

## Recombinant Human IL-16 Protein, HEK293 Lysate

**Catalog Number:** 230-10272-100

### Data Sheet / Certificate of Analysis (CoA)

Last Revised: 11/26/2024

### Source

---

<b>Species</b>	Human
<b>Accession Number</b>	Q14005-2
<b>Gene Symbols</b>	IL16
<b>Synonyms</b>	Interleukin-16, transfected HEK293 cell lysate, Interleukin 16, Lymphocyte Chemoattractant Factor, Prointerleukin 16, Interleukin 16 (Lymphocyte Chemoattractant Factor), Neuronal Interleukin 16, Pro-Interleukin-16, PrIL-16, PRIL16, NIL16, LCF.
<b>Expressed Region</b>	Pro1204-Ser1331

### Preparation

---

<b>Expression System</b>	HEK293 cells
<b>Tag</b>	N-terminal histidine tag
<b>Purity</b>	Unpurified cell lysate. HEK293 cells transfected with expression vectors harboring target gene were harvested and washed with PBS twice. The cell pastes were re-suspended with ice-cold PBS containing mammalian cell protease inhibitor cocktail and further lysed with freeze-thaw cycles. After clarifying with 20,000 g centrifugation at 4°C for 30 min, the lysate was aliquoted, lyophilized, and stored at -80°C immediately. Protein concentration was determined by BCA kit (Thermo Scientific, Inc.) using BSA as protein standard. The gene overexpression in lysate was confirmed by Western blotting using anti-His tag antibody and/or target-specific antibodies and the lysate derived from HEK293 cells transfected with the empty expression vector was used as a negative control.
<b>Molecular Weight</b>	Recombinant protein has a calculated molecular mass of 13. The actual molecular weight may increase differently due to the potential post-modifications (PTMs).

### Specifications

---

<b>Format</b>	Lyophilized powder
---------------	--------------------

<b>Formulation</b>	Lyophilized from a 0.2 $\mu\text{m}$ filtered solution in PBS containing mammalian cell protease inhibitor cocktail
<b>Concentration</b>	Determined by Pierce BCA protein assay kit
<b>Reconstitution</b>	Briefly spin the vial and bring the contents to the bottom prior to opening. It is recommended to reconstitute at 0.5 - 1 mg/mL with sterile deionized water.

## Shipping

---

The product is shipped with ice packs. Upon arrival, immediately store it at the temperature recommended below.

## Storage/Stability

---

- Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
- Upon arrival, the lyophilized protein may be stored for 2 weeks at 4°C.
- For long term storage, it is recommended to store desiccated below -20°C in a manual defrost freezer.
- Generally, the shelf life is up to 12 months from date of receipt at -20°C or -80°C under sterile conditions.
- Following reconstitution, the protein may be stored for 2 weeks under sterile conditions at -20°C. For long term storage, it is recommended to make appropriate aliquots at -20°C or -80°C.