
Recombinant Human Ubiquitin B Protein, HEK293

Catalog Number: 230-30221

Data Sheet / Certificate of Analysis (CoA)

Last Revised: 11/27/2024

Source

Species	Human
Accession Number	P0CG47
Gene Symbols	UBB
Synonyms	Polyubiquitin-B, Ubiquitin, UBB
Expressed Region	Gln2-Gly228

Preparation


Expression System	HEK293 cells
Tag	N-terminal histidine tag
Purity	>90%
Molecular Weight	Recombinant protein product has a calculated molecular mass of 26 kDa. The partial purified product cleaved into 12 and 14 kDa two fragments.
Endotoxin Level	<0.5 EU per ug of the protein as determined by the LAL method

Specifications

Format	Lyophilized powder
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS (pH 7.4) with 1% mannitol and 5% trehalose
Concentration	Determined by Pierce BCA protein assay
Reconstitution	Briefly spin the vial and bring the contents to the bottom prior to opening. It is recommended to reconstitute at 0.5 - 1.0 mg/mL with sterile deionized water.

SDS-PAGE Image

Figure 1. Deglycosylation analysis of purified recombinant proteins. Lane 1, protein ladder (kDa). Lane 2, glycosylated protein (delivered to customers). Lanes 3 and 4, deglycosylated proteins under native and denature conditions respectively.



Shipping

The product is shipped with ice packs. Upon arrival, immediately store it at the temperature recommended below.

Storage/Stability

- Use a manual defrost freezer and avoid repeated freeze-thaw cycles.
- Upon arrival, the lyophilized protein may be stored for 2 weeks at 4°C.
- For long term storage, it is recommended to store desiccated below -20°C in a manual defrost freezer.
- Generally, the shelf life is up to 12 months from date of receipt at -20°C or -80°C under sterile conditions.
- Following reconstitution, the protein may be stored for 2 weeks under sterile conditions at -20°C. For long term storage, it is recommended to make appropriate aliquots at -20°C or -80°C.